

Official site

[colocolo-brush.com](http://colocolo-brush.com)

Patented Manufacturer T.W.C Inc.  
1-3-17, Katayama-cho,  
Suita-shi, Osaka, Japan  
+81-66368-7117

Surprisingly

effective

roller brush

COLOCOLO  
BRUSH



# One-Stop Oral Care

Achieve everything you need for your oral care.

A single Colocolo brush will make your life positive.  
Please feel improvement on your teeth, mouth, mind and body.

## ONE-STOP ORAL CARE

Gum  
Disease

Beauty

Health

Provide new-style tooth brushing: place the brush to your teeth and roll along them. This brushing style will not only prevent and cure gingivitis, pyorrhea, and bad breath, but also make you feel relaxed after massaging the gums with the brush and improve your beauty including weight loss. Furthermore, strengthening the immune system and stabilizing blood pressure by using this brush are scientifically proven. Why don't you start one-stop oral care that supports your healthy life with a Colocolo brush?

This brush  
registered  
by Japanese  
oral dentist  
association

The effectiveness  
of critical test  
and study was  
taken by  
the association.


# Gum Disease

- Gingivitis
- Pyorrhea
- Bad breath

## The idea of rolling a toothbrush can keep your gums healthy.

Rolling a Colocolo brush inside the mouth will stimulate your gums and increase their blood flow like a massage. As a result, tighten the gums around the teeth, keep them healthy, and prevent and cure gum disease.




The case study proves that patient gum  part swore.

Used Colocolo bush and after 3 days, patients same gum area shrunken and after 6 months entire gum area firm and rejuvenated.

## You experiment improvement of gum area in 7 days.

This brush gently massages on the gum area resulted in avoid scratch or damaged gum. Within the just one week, gum gets rejuvenate and firm feeling experienced by the user.




\*The illustration is an image

## Stabilizing your oral bacteria level helps prevent and cure your bad breath.

A coated tongue appears on the upper surface of your tongue like a white moss. Although seeing the slightly white surface of your tongue is harmless, a thick white coating causes bad breath. Rolling a Colocolo brush along your teeth and gums helps stabilize your oral bacteria level and naturally reduce this white coating. You can roll the brush along the upper surface of your tongue too. Please use the Colocolo brush for preventing and improving your bad breath.



Tung area seems improved by using this brush.



# Beauty


- Whitening
- Healing
- Relaxation
- Diet



## Restore Your Original White Pearls for a Healthy Smile

The Colocolo brush vitalizes all oral health in result, bring your teeth healthy white. The dental and oral health association recognized the effectiveness of Colocolo brush. Removes tobacco tea and coffee stain, whiten your teeth.

Relaxing effect by brushing




## Pleasant effects on your mind and body

The data shows that massaging the gums with a Colocolo brush increases alpha brain waves and stabilizes the blood pressure and mood. Many beauty benefits, such as healing, relaxing, and weight loss, can be expected by using this brush.

# Health

- Immune improvement
- Stable action of blood pressure


Change of blood pressure by Colocolo brush



## Building healthy blood

Using a Colocolo brush helps increase blood flow to the gums. That can be expected to control your blood sugar and blood pressure. Significant improvements in blood health were seen after those with chronic hypertension used this brush.

Relationship between salivary volume and oral bacteria count



## Increasing saliva production

Brushing your gums and massaging the inside of your cheek with the outer side of the brush head stimulates salivary glands and increases saliva production. As saliva contains antibacterial compounds, this will help prevent tooth decay, gum disease, and bad breath.

## FUNCTIONALITY

# A Colocolo brush can fit into the teeth about 10 times better than conventional tooth brushes

Brushing regularly with the Colocolo brush can help solve your oral problems.



### POINT.1

## A 360-degree Roller Head

It is equipped with the world's first "personal Spin Disc Technology." Each disc brush has a diameter of 17.5 mm and about 650 nylon strands welded ultrasonically, radiating out from the center. The circular brush made of 14 discs, put together using a patented method, rolls over the teeth, peeling, absorbing, and removing plaque.

### POINT.2

## Over 9,000 Tips of Hair

There are over 9,000 with just the right softness. It is about 10 times that of regular toothbrushes and can brush between the teeth without forcing its way through, proving the gums with just the right amount of simulation.

### POINT.3

## Easy-grip makes a good result

It has an easy-to-hold grip with a universal design based on human factors and ergonomics. Hold the grip like you would a pencil and roll the brush over without applying pressure. Doing just that will give you a significant improvement.

Presenting clinical results at the dental society



evidence is here

### POINT.4

## Lightly Roll It Around

Place it on your teeth and roll it 30 round-trips back and forth. Lightly rolling it around removes about 90% of plaque. It is a total oral care solution, taking care of the surfaces of your teeth, spaces between your teeth, gums, and the top of your tongue.

### POINT.5

## Acquired Patent

The ColoColo Brush is the only toothbrush that has acquired a basic patent in the manufacturing method and device for producing 360-degree brushes.

Patent No. 3981290  
Patent No. 4037653  
Patent No. 4037739


\* There are products with similar concept and shape from other brands. However, they have no relations to our Colo Colo Brush whatsoever. Please make sure to check when purchasing.

USAGE & INDICATION

# How to Use a Colocolo Brush & its Beneficial Effects




## Rolling Around & Removing Plaque

Rotating brush creates ultrasonic vibration in result removes hard plaque from your teeth without damaging your teeth or gums. Roller action superior than a regular toothbrush.



## Moving a Colocolo brush along the teeth like using a spatula can powerfully remove plaque.



By splitting the brush head right after placing the brush to your teeth, the brush can easily reach through the gaps between your teeth and periodontal pockets. As the brush rolls, its head rotates more than 9,000 extra fine bristles. That can powerfully remove plaque like rubbing the surface of the teeth with a spatula.



Just placing the brush to your teeth and rolling up and down 30 times can remove an overwhelming amount of plaque, which is not possible with conventional tooth brushes.

## Massage Effect on your gum.

Stimulate and massage gums at points rather than faces to promote blood flow.



## Squeaky Clean & Feels Great! For Surprisingly Smooth Teeth

Touch the surfaces of your teeth after brushing. They should be squeaky clean. Most Colo Colo Brush users have experienced their teeth becoming smoother.

Subjects: A total of 35 men, 18 males and 17 females aged 9-12. Experimental method: The upper chin 4 front teeth which 100% plaque adhesion was recognized at the time of plaque staining were polished sideways with a roller brush, and the plaque residue was checked for each of the 10th and 30th reciprocations.

I recommend to use a Colocolo brush.

To keep a smile on my face,  
I cannot live without  
this brush.

Actress **Miho Matsubara**

As an actress and a singer, I believe that it is important for me to make my audience happy and for that, I must feel good. That is why I pay attention to beauty and health. One of the things that I always rely on is the Colo Colo Brush. All I have to do is to roll it over my teeth and gums to make my teeth smoother and make my gums feel better. I also feel that my gums tighten around my

teeth again. I feel more confident acting on film or stage when my teeth look good. What makes me happier than anything is that I can see my audience with a big smile. I wouldn't be able to start my day without it for my beauty and health! I would like to keep myself up and running alongside the Colo Colo Brush so that I can smile and meet more people.



#### PROFILE

Born in the year 1976. Debuted Japanese idle grope "Kuno Ichi". Recorded title Love Gipsy/Like Some with CD. Acted in many Japanese TV drama and movies in the past. Recent single recording "Dream of Miho. Talent skill with Judo. Japanese folkloric Dance and Jogging is my hobby.

## Even Dentist open eyes!

Why do people get cavities even when they brush their teeth every day? We believe that the problem is the toothbrushes. We needed a new toothbrush that could do what regular toothbrushes couldn't do. That is how we Colo Colo Brushes came to be. It removes a significant amount of plaque, has brush tips that reach into the dental pockets, massages the gums for better blood flow, and brushes the tongue. Our one-stop oral care device has been clinically tested for effectiveness and has received great reviews from those that suffer from many oral health conditions.



Patented

Patent No. 3981290

Color Variation



BLUE

PINK

GREEN

YELLOW

GREY

**MADE IN JAPAN**

Product Name	ColoColo Brush
Model	Ver. 2017
Brush specifications	over 9000 nylon hair
Color	Blue Pink Yellow Green Gray
Material	Handle: ABS resin (heat-resistant 70 ° C) Hair: Nylon (made by Dupon) 0,1 mm
body size	H17,5cm W1,4cm D2,5cm
Product weight	12,5g
Country of manufacture	Japan
Patent number	No. 3981290
Design number	No. 1549462

I recommend to use a  
Colocolo brush.



Bio Dental Clinic Representative  
Ohara Dentistry Clinic Director  
Aichi Gakuin University,  
Matsumoto Dental University,  
Periodontal department  
Part-time lecturer

Doctor of dentistry  
**Ohara Morikatsu**

Simply brushing your teeth is not enough for proper oral care. "Keep your mouth clean and healthy". It is very important to be aware of this when you take care of your teeth and mouth. This Colocolo brush is innovatively designed so that anyone can easily perform proper oral care. Nothing but gently stroking every area of your mouth, including your teeth, gums, tongue, and the inside of your cheek, will help stimulate saliva flow and keep your mouth clean and healthy. Achieve a healthy oral environment by changing your approach and your tool for oral care. Keeping your mouth clean and healthy will make you feel happy more and more every day.